

Hephaestus V V094, V213, & V900 Vacuum Drying & Purging Ovens

The Hephaestus vacuum oven is designed for drying and vacuum purging materials which are thermosensitive, oxidative or easily decomposable. This includes cannabis flower just as much as shatter or isolates.

To further reduce oxidation e.g. to protect THCa you can purge or vent the oven with inert gas to displace any oxygen present.






Safety:

- ETL and CE rated
- Overtemperature Alarm

Features & Benefits:

- Benchmark heat up time – 40% faster
- 304 stainless steel with mirror polished surface for easy cleaning and maintenance
- A PID controller with over temperature alarm and timing function ensures precise and reliable control
- Large dual layered tempered glass door with LED lighting providing safety and good observation of your drying process
- Smart door adjustment system with silicone door gasket for good vacuum performance
- Inert gas vent valve
- LED lights
- KF25 Vacuum Port and Barbed Purge Port



Hephaestus V V094

HGS SKU Number	BLU-VOV-DZF6094-US
Power Specification	120V 60Hz 1ph – 1350W
Chamber Size W×D×H (mm)	90L (450x450x450)
Temperature Range	ambient + 10°C to 200°C
Shelves	4
Dimensions (cm)	88 x 76 x 91
G.W (KG)	166





Hephaestus V V213

HGS SKU Number	BLU-VOV-DZF6213-US
Power Specification	120V 60Hz 1ph – 1550W
Chamber Size W×D×H (mm)	216L (600x600x600)
Temperature Range	ambient + 10°C to 200°C
Shelves	3
Dimensions (cm)	106 x 97 x 105
G.W (KG)	265





Hephaestus V V900

HGS SKU Number	BLU-VOV-DZF6930-US
Power Specification	230V 60Hz 1ph – 5800W
Chamber Size W×D×H (mm)	900L (750x1160x1050)
Temperature Range	ambient + 10°C to 200°C
Shelves	5
Key Technical Specs	Display Resolution: 0.1 °C Temperature Stability: ±1°C Vacuum: 1.3 mbar
Dimensions (cm)	222 x 165 x 130
G.W (KG)	598