Benchmark E



Prefilled Tubes and Bulk Beads: The Right Bead for Every Sample

- Wide range of kits available
 Variety of bead sizes and types to
 meet your sample needs
- __ TriplePure™ Quality
 Nuclease, Protease and Nucleic
 Acid free
- Prefilled 2.0ml tubes or bulk beads
 Convenience of prefilled tubes
 Bulk beads for filling your own

Benchmark Scientific's homogenizer bead kits come in a variety of materials and sizes for processing a wide range of samples. Prefilled for convenience, the unskirted 2.0ml tube kits will fit in all popular bead homogenizers including Benchmark Scientific's BeadBug™ and BeadBalaster™ 24. Tested to ensure they are RNase, DNase, protease, and nucleic acid free, the TriplePure Zirconium Bead Kits are supplied in five bead sizes to homogenize everything from bacteria to large tissues. Acid washed silica and stainless steel bead kits as well as specialty bead mixes and bulk beads are also available.

Bead material	Bead Diameter	Prefilled 2.0ml tubes, pk/50	Bulk cat. no./bottle		
TriplePure Zirconium	0.1mm	D1032-01	D1132-01TP/250g		
TriplePure Zirconium	0.5mm	D1032-05	D1132-05TP/250g		
TriplePure Zirconium	1.0mm	D1032-10	D1132-10TP/250g		
TriplePure Zirconium	1.5mm	D1032-15	D1132-15TP/250g		
TriplePure Zirconium	3.0mm	D1032-30	D1132-30TP/300g		
TriplePure Zirconium	0.1, 0.5, 1.0, 1.5, 3.0mm	D1032-SK*	NA		
Silica	0.1mm	D1031-01	D1131-01/200g		
Silica	0.5mm	D1031-05	D1131-05/200g		
Silica	1.0mm	D1031-10	D1131-10/200g		
Stainless Steel	2.8mm	D1033-28	D1133-28/1,000 balls		
Garnet, Zirconium	Garnet shards/6mm zirconium satellite	D1033-30G	NA		
Silica, Zirconium	Mixed size	D1034-MX	NA		
Empty tubes, with screv	D1031-T21				
Empty tubes, with screv	D1031-T20				
BeadBug Bead Homoge	D1030				
BeadBlaster 24 Bead Ho	D2400				
*Starter kit - contains 10 of	each size.				

	Bead material	Silica			Zirconium				Steel	Garnet with Zr Satellite	Silica and Zirconium	
	Bead diameter	0.1mm	0.5mm	1.0mm	0.1mm	0.5mm	1.0mm	1.5mm	3.0mm	2.8mm	0.5/6mm	Mixed
Animal	Adipose							Х				
	Artery							Х	Х			
	Bone									Х		
	Brain							Х	Х			
	Insect - fly								X	Х		
	Insect - Tick								X	X		
	Insect - Mosquito						Х	Х	X			
	Gonad							X	X			
	Hair								X	Х	Х	
	Heart							Х	X			
	Kidney							Х				
	Liver							X	X			
	Lung							Х	X			
	Muscle								X	Х		
	Pancreas							Х	X			
	Skin								X		Х	
	Spleen							Х	X			
	Teeth								Х	Х		
Microbial	Algae		Х			Х						
	Bacteria	Х	X		X	Х						
	Fungi		Х	Х		Х	Х					
	Mold		Х			Х						
	Spores		Х	Х		Х	Х					
	Yeast		Х	Х		Х	Х					
Plant and Soil	Leaves								Х	Х		
	Nuts									Х		
	Plant, soft tissue								Х	Х		
	Plant, tough tissue								Х	Х	Х	
	Seeds								Х	Х		
	Stems								Х	Х		
	Soil							Х				Х
	Coral									Х		
Fecal	Stool							Х				Х
	Scat								Х	Х		Х

Above is a partial list of sample types and suggested beads. For additional information and technical support, contact Benchmark Scientific at info@BenchmarkScientific.com.



PO Box 709 Edison, NJ 08818 TEL: 908 769-5555 Fax: 908 222-1864 email: info@BenchmarkScietific.com Web: www.BenchmarkScientific.com