



➤ Technical data

Hei-TORQUE Core

Supply voltage	100 - 240 V (50/60 Hz)
Rated power (W)	105
Protection class (IEC 61140)	1
Protection class (IEC 60529)	IP42
Sound pressure level (dB(A)) (based on IEC 61010)	< 70
Drive	EC motor
Motor protection	Software monitoring with error messages
Overload	Auto switch-off
Rated torque (Ncm)	40
Max. torque (Ncm)	65
Performance display	Yes
Speed range (rpm)	20 - 2,000
Speed setting	Rotary knob
Speed accuracy (%)	± 1
Speed control	Electronic
Permissible duty cycle	Continuous operation
Min./max. impeller shaft Ø (mm)	0.5 – 10.5
Chuck	Keyless
Interfaces	No
Timer	Yes
Timer range (min)	0 - 6001
Timer accuracy (%)	± 1
Control panel	TFT display, monochrome, 2"
Max. volume H ₂ O (L)	25
Max. viscosity (mPa s)	10,000
Number of gear steps	1
Dimensions (mm)	W 70 H 281.5 D 195
Weight (kg)	2.3
Support rod (Ø x L) (mm)	13 x 160
Weight of support rod (kg)	0.165