




# Polycarbonate Serum Vials

## DATA SHEET



Vial Dimensions	Catalog #	RPC0010	RPC0005
	Size / Volume	10 ml	5 ml
	A • Diameter	23.7 mm	22.3 mm
	B • Height (w/o Cap)	62 mm	52 mm
	C • Height (w/ Cap)	64.1 m	54.1 mm



TriForest Serum Vials are ideal for sampling, storage, production and packaging of steam sterilized biologicals. They are ideal for leak-proof shipment of reagents and buffer solutions. All vials and caps have lot number traceability. These shatterproof, resistant containers are useful for packaging serum, veterinary, pharmaceuticals, cosmeceutical, and nutraceutical extracts. They can be used with most aromatic compounds and essential oils. They are provided non sterile with caps separately packed.

TriForest Enterprises, Inc.

6 Chrysler, Irvine, CA 92618, USA | (888) 777-9956 | (949) 380-8866 | Fax: (949) 380-9955 | [www.triforest.com](http://www.triforest.com) | [info@triforest.com](mailto:info@triforest.com)  
© Copyright of TriForest Enterprises, Inc. Vol. 1

# Polycarbonate Serum Vials

## VIAL & CAP SPECIFICATIONS

Vial Specification	Catalog #	RPC0010	RPC0005
	Overhead Space	20%	
	Transparent	YES	
	Leakproof	YES	
	Vials Shrink Wrapped in Tray	YES	
	Particle Free	NO	
	Sterile (Gamma Irradiation)	NO	
	Temperature Range	-125° C to +125° C	

Sterilization	Autoclavable	YES	
	Maximum Autoclave Temperature	125° C	
	Autoclave Cycles	≤ 10	

Cap Specification	Size	20-415	20-415
	Cap Type	Standard Seal Cap	
	Standard Color	Blue	
	Material	Polypropylene	

Packaging	Number of Vials / Pack	100	100
	Number of Trays / Case	16	16
	Total Number of Products / Case	1600	1600
	Case Dimension (L x W x H inches)	19.5 x 19.5 x 14	18.5 x 18.5 x 12.5
	Case Dimension (L x W x H cm)	50 x 50 x 36	47 x 47 x 32
	Case Gross Weight	40 lbs / 17.5 kg	33 lbs / 15 kg