

labForce

powered by *Thomas Scientific*



TS-B Laboratory Balances

BASIC PORTABLE WEIGHING

TS-B Series

Laboratory Balances

Five-Year Warranty. Worry-free for five years, all replacement costs and labor included.

Portable. Operate in the field or move from lab-to-lab. All TS-B Series come complete with 12v DC Power Adapter or insert six AA (batteries not included) for portable applications. Automatic shut off mode conserves battery life and the low battery indicator ensures sample weighing is not interrupted.

Simple Operation. Four dedicated keys ON/OFF, CAL, UNIT and TARE are ideal for student use.

Oversized Draft Shield with a stainless steel lift off cover is standard on the TS-B-123. The internal dimensions (150 x 140 x 107mm) provides easy access and ample space when working with larger samples. The TS-B-123 also includes a 100g calibration weight.

Plastic Weigh Bowl. The TS-B Series (except for the TS-B-123) includes an oversized weigh bowl to allow for easy weighing of items.

Stainless Steel Weigh Pans are easy to remove for cleaning.

Rugged, Built to last. Constructed to endure environments from teaching labs to light industrial use.

Priced Right. The LabForce TS-B Series Balances are the solution for organizations with limited budgets.

Five year warranty included on all models

Model	TS-B-123
Maximum capacity	120 g
Readability [d]	0.001 g
Taring Range	0-120 g
Repeatability (s)	0.001 g
Linearity	0.003 g
Adjustment	External
Stabilization Time	3 sec.
Pan Size	70mm (7.76in)
Unit Dimensions	225 x 177 x 157mm (8.9 x 6.7 x 6.2in)
Catalog Number	1160U14



Model	TS-B-502	TS-B-1502	TS-B-3001	TS-B-6000
Maximum Capacity	500 g	1500 g	3000 g	6000 g
Readability [d]	0.01 g	0.05 g	0.1 g	1 g
Taring Range	0-500 g	0-1500 g	0-3000 g	0-6000 g
Repeatability	0.01 g	0.05g	0.1g	1g
Linearity	0.02 g	0.1 g	0.2 g	2 g
Adjustment	External	External	External	External
Stabilization Time	3 sec.	3 sec.	3 sec.	3 sec.
Pan Size	145 x 145mm (5.7 x 5.7in)	145 x 145mm (5.7 x 5.7in)	145 x 145mm (5.7 x 5.7in)	145 x 145mm (5.7 x 5.7in)
Unit Dimensions	225 x 177 x 70mm (8.9 x 6.7 x 2.8in)	225 x 177 x 70mm (8.9 x 6.7 x 2.8in)	225 x 177 x 70mm (8.9 x 6.7 x 2.8in)	225 x 177 x 70mm (8.9 x 6.7 x 2.8in)
Catalog Number	1160U15	1160U16	1160U17	1160U18

