

Sonifier Cell Disruptor



SFX 150 Microtips

Microtips for the SFX Sonifier are available in four sizes to handle volumes up to 150 ml. The microtips are also designed to process samples in the smallest diameter vessels.

The 1/8" microtip is already provided with the SFX Sonifiers.

Part #	Processing Volume	Tip Diameter	Approx. Length (working)
109-122-1066	1ml & under	3/32" (2.4mm)	4"
109-122-1065	3 -10 ml	1/8" (3.2mm)	4"
109-122-1182	5 -25 ml	3/16" (4.8mm)	4"
109-122-1064	10 - 150 ml	1/4" (6.4mm)	4"
109-116-1566	3-10 ml (4) tips	1/8" (3.2mm)	4"

Note: The 109-116-1566 has (4) 1/8" microtips on a 1" center horn.

Note: These tips are for use with the Branson 150 watt SFX and SLP Sonifier only.

Sonifier Cell Disruptor

Global Technologies. Local Solutions.

Branson's unmatched global resources ensure optimal solutions for the most challenging ultrasonic precision cleaning and materials joining applications. With 1600+ employees in over 70 sites worldwide, we can rapidly respond to our customers' needs, wherever they are located.

Branson understands local markets and regulations, and the open collaboration among our global offices and extensive staff of application specialists quickly produces solutions for customers. We partner with companies of every size and

scope to help resolve critical issues ranging from market changes, product quality/life cycles and production costs, to employee safety and environmental compliance. And, Branson's commitment to the publicly funded industrial collective research initiative and cooperation with various research centers have produced widespread innovations in ultrasonic cleaning and plastics joining technology that are benefiting companies everywhere.

Branson is a leading innovator in Emerson, a diversified global

manufacturing and technology company. Emerson is a progressive company working to develop products and processes to resolve global issues, including energy supply and distribution, information and networking, climate and environment, manufacturing efficiency, and product performance. Emerson is dedicated to the pursuit of breakthrough technology developments never before envisioned.

Branson brings exceptional quality, fast delivery and competitive prices. Contact your nearest regional center today.

BRANSON

Branson Ultrasonics Corporation
41 Eagle Road, Danbury, CT 06813-1961
(203) 796-0400 • Fax: (203) 796-9838
www.bransonultrasonics.com
e-mail: info@bransonultrasonics.com

SFX Microtips © Branson Ultrasonics Corporation 2016
The Emerson logo is a trademark and service mark of Emerson Electric Co.
Revised and printed in the U.S.A. 1/2016