

The **NEW** Transferpette® S

A perfect fit for every hand!

BRAND. For lab. For life.®



Precision workmanship, high grade materials, "Made in Germany" quality – it's what you see on first glance and appreciate every time you use the new Transferpette® S. The finer points of the design will be evident when the instrument is in your hand. Experience how simple, easy and efficient pipetting can be with the new Transferpette® S. Single- or multi-channel, there's a Transfepette® S pipette for your application.

- **Slim hand grip**—Provides a comfortable and secure fit in any hand
- **One-handed operation**—Perform all functions (set volume, pipette, tip ejection) without changing your grip
- **Volume-change protection**—Simple to operate and clearly shows when the volume setting is locked
- **Slim shaft**—Narrow tubes? With the slim shaft, easily pipette into narrow tubes such as 15 mL conical tubes
- **Easy cleaning**—The shaft tucks inside the handgrip for smooth surfaces, making cleaning easier
- **Low forces**—Short pipetting stroke, easy mounting and tip ejection make pipetting a breeze
- **Universal tip cone** fits tips from most quality manufacturers
- **Completely autoclavable** at 121 °C
- **Color-coded pipette** accents make pipette and tip selection easy
- **Easy mounting and ejection** of tips with O-rings on multi-channel pipettes
- **3-year warranty** for single-channel pipettes



Reaktionsgefäßständer
Micro Test Tube Rack
BRAND

Thomas
Scientific

BRANDTECH
SCIENTIFIC, INC.

Transferpette® S

Additional product information can be found at ThomasSci.com



Transferpette® S pipettes, adjustable volume (fixed volume also available)

	Capacity, µl (color-coded)	Description	BrandTech® No.	Thomas No.
●	0.1–1	D-1	705868	20A00F974
●	0.1–2.5	D-2.5	705869	20A00F983
●	0.5–10	D-10	705870	20A00F982
●	2–20	D-20	705872	20A00F981
●	5–50	D-50	705873	20A00F980
●	10–100	D-100	705874	20A00F979
●	20–200	D-200	705878	20A00F978
●	100–1000	D-1000	705880	20A00F977
●	500–5000	D-5000	705882	20A00F976
●	1000–10000	D-10000	705884	20A00F975

Transferpette® S-8 and -12 multi-channel pipettes

	Capacity, µl (color-coded)		BrandTech® No. 8-channel	Thomas No. 8-channel	BrandTech® No. 12-channel	Thomas No. 12-channel
●	0.5–10	M8-10	705900	20A00F990	705920	20A00F995
●	5–50	M8-50	705906	20A00F994	705926	20A00F999
●	10–100	M8-100	705908	20A00F993	705928	20A00F998
●	20–200	M8-200	705910	20A00F992	705930	20A00F997
●	30–300	M8-300	705912	20A00F991	705932	20A00F996



Pipette sets and packages, Transferpette® S, adjustable volume, includes tips and shelf/rack (w/sets), benchtop rack (w/packages)

Type	Description	BrandTech® No.	Thomas No.
Micro	D-2.5, D-10, D-100	705890	20A00F984
Mid	D-20, D-200, D-1000	705891	20A00F989
Macro	D-1000, D-5000, D-10000	705892	20A00F988
Standard	D-10, D-100, D-1000	705893	20A00F987
Pipette package 1	D-2.5, D-10, D-100, D-200, D-1000	705896	20A00F986
Pipette package 2	D-10, D-100, D-200, D-1000, D-5000	705897	20A00F985

