



4 Research Court, Suite 300
Rockville, MD 20850
877-411-2041
Services@ibtbioservices.com

Human anti-*C. difficile* toxin A

Catalog #: 0250-001

Lot #: 1303004

Immunogen: Recombinant *C. difficile* toxin A.

Description: A human IgG produced in *N. benthamiana* and is reactive to *C. difficile* toxin A.

Supplied: 100 µg is supplied in PBS pH 7.4 at a concentration of 2.09 mg/mL.

Purification: Antibody is purified using immobilized protein A.

Clonality: Human variable, human constant of the IgG₁ isotype.

Relevance: This monoclonal is a potent neutralizer of toxin (IC₅₀ = 1.0 ng/ml) *in vitro*.

Recommended Dilutions:

ELISA: Recommended antibody concentration from 10 – 100 ng/mL.

WB: This antibody has not been tested in Western Blot analysis.

Storage: 2-3 weeks +4°C, -80°C long term. Avoid repeated freeze/thaws.

Cross Reactivity: No known cross reactivity.

Related Products:

Human anti- *C. difficile* toxin B Antibody

Catalog # 0250-002

Anti- *C. perfringens* ETX Antibody

Catalog # 0240-001

Intended for research use only, not for human, therapeutic, or diagnostic applications.