

Product Name:	LF171WWW/0HC
Description:	<p>General Purpose Auto Defrost Hydrocarbon Laboratory Freezer</p> <p>17 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Mechanical thermostat</p> <p>Adjustable operating temperature range: -15°-25°C</p> <p>Three adjustable shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>One year parts and labor warranty; plus additional four year compressor parts warranty</p> <p>Exterior dimensions: 28" W x 29 3/4" D x 68"H</p> <p>Auto defrost</p> <p>115V, 60Hz, 5 Amps, 1/4 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 227 lbs.</p> <p>DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p>