



QP-FORM-1650-160-#9800
Ver#4 Prepared by GA approved by MAF6/19/13
F60-WIDRG-000140-VPS (9800) Specification Sheet

SPECIFICATION SHEET
Item# F60-WIDRG-000140-VPS
Short Name # 9800

ITEM DESCRIPTION:

The 0-14 Hydrion Indicator Strips consists of 100 pH strips with a color chart with color matches at each full pH increment from 0-14. (0-1-2-3-4-5-6-7-8-9-10-11-12-13-14)

UNIT

Each unit shall contain one flip-top vial with color chart label and 100 plastic pH strips.

ITEM SPECIFICATIONS

- (1) **Vial** – White plastic vial with attached flip top.
- (2) **COLOR CHART** – With color matches at 0-1-2-3-4-5-6-7-8-9-10-11-12-13-14.
- (3) **pH Test Strips** – The vial shall contain 100 pH test strips, each measuring 1/4” x 2.75”

DIRECTIONS

Remove plastic strip from vial and dip into the solution being tested. Strip should change color instantly. Remove at once and compare the resulting color to the matching pH color chart.

