



Intelli Mixer RM-2M, Medium



Product Description

Powered by V-Spin – the world’s most advanced, patented mixing and liquid-handling technology – the Intelli-Mixer RM-2 is truly the latest in cutting-edge science. V-spin harnesses the energy of an unequaled stepper motor, allowing for a wide variety of mixing, vibrating, and rotating cycles. This multipurpose unit is a rotator, vortex mixer, and blotter all in one. Designed to be cold room safe, and with an unprecedented 34 modes of operation, this versatile piece is unrivaled in today’s market.


Features

- Rocker, Vortex Mixer, and Blotter
- 34 automated intelligent mixing programs
- 3 Year Warranty

Specifications

Rotation Speed	1 - 99 RPM
Rotation Speed Settings Step	1 RPM
Overall Number of Programs	28 Standard + 3 Custom
Environmental Temperature	10 - 45 °C
Relative Air Humidity (at 20 °C)	up to 80 %
Power Supply Adapter	DC 12V, 1.5 A

Accessories

	IMIX-16 Standard Racks Tube Size: 9 mm Number of Tubes: 30		IMIX-15 Standard Racks Tube Size: 11 mm Number of Tubes: 32
	IMIX-14 Standard Racks Tube Size: 13 mm Number of Tubes: 28		IMIX-13 Standard Racks Tube Size: 16 mm Number of Tubes: 16
	IMIX-12 Standard Racks Tube Size: 30 mm Number of Tubes: 10		IMIX-17 Mix Racks # of 11mm Tubes: 8 # of 13mm Tubes: 8 # of 16mm Tubes: 4 # of 30mm Tubes: 2
	IMIX-18 Rocking Platform		IMIX-19 Blotting Rack Tube Size: 50 ml Holds 6 Tubes

