

## Vaccinia NYCBOH Antibody

**Catalog #:** 01-0003

**Immunogen:** Vaccinia (NYCBOH strain) infected cell lysate supernatant

**Description:** Rabbit polyclonal for Vaccinia vaccine

**Raised in:** New Zealand White Rabbits

**Purification:** Neat serum

**Clonality:** Polyclonal

**Buffer:** serum + 0.01% Azide

**Size:** 250ul

**Relevance:**

Vaccinia represents a surrogate strain for the Category A (NIAID) and HHS Select Agent variola major.

**Recommended Dilutions:**

ELISA: Assay-dependent dilution


WB: Assay-dependent dilution. If probing against infected cell lysate, cross-adsorption recommended.

IF: No Data

**Storage:** 2-3 weeks +4°C, -20°C long term

**Species Reactivity:**

Does not react with Venezuelan Equine Encephalitis TC-83, Chikungunya 181/25, Junin candid #1, or Rift Valley Fever MP12. Some cross-reactivity observed with fetal bovine serum and BSC-40 cells.



-Serum tested against infected cell lysate and cells alone (mock)  
-Serum diluted 1:1,000 and incubated with 2% FBS for 2hr  
Prior to blotting

**Intended for research use only, not for human, therapeutic, or diagnostic applications.**