

NXF3 Polyclonal Antibody

Catalog number: 13275-1-AP


Size: 69 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

NXF3; nuclear RNA export factor 3, NXF3, TAP like protein 3, TAPL 3, TAPL3



Immunohistochemical of paraffin-embedded human gliomas using 13275-1-AP (NXF3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using 13275-1-AP (NXF3 antibody) at dilution of 1:50 (under 40x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human; other species not tested.
Calculated NXF3 MW:	531aa,60 kDa
Observed NXF3 MW:	
Positive IHC detected in	Human gliomas tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3998
GenBank accession number:	BC031616
Gene ID (NCBI):	56000
Full name:	Nuclear RNA export factor 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.