

Product Name:	ABS VS 3
Description:	<p>ABS VS Class Vapor Shipper Includes 1 Canister Total Vial Capacity (6/Cane): 102 Liquid Nitrogen Capacity: 3.6 Liters ½ CC Straws (1 Level Bulk): 440 ½ CC straws (10/Cane): 280 Static Evaporation Rate (Liters/day): 0.26 Static Holding Time (days): 13 Neck Diameter (mm): 69.9 Exterior Height (mm): 467.4 Exterior Diameter (mm): 221 Canister Height (mm): 279.4 Canister Diameter (mm): 66.6 Optional protective shipping container ensure upright shipping High strength neck tube reduces liquid nitrogen loss Lid with locking provision Hydrophobic liquid nitrogen absorbent system 2-hour charging capability One-year parts and labor warranty One-year vacuum warranty</p>