



Note: labeling is forbidden on the Freezing Bags body;  
a small label can be placed on the surface of the inner PE bag  
and a product label on the outer PE bag

LePure Biotech Co., Ltd. Copyright. All rights reserved.

Item No.	Qty	Description	Inventory Code	Material
D1	1	500mL Freezing Bags(0.2)	ACAX0022	EVA
A1	2	0.203"ID/0.280"OD(5.15mm/7.1mm)EVA	EABX0083	EVA
B1	1	2.95*4.15mm EVA Tubing,Length 18cm	DDAX0063	EVA
C1	1	Freeze tube tubing assembly	DDAX0060	PVC/Silicone/PE
F1	1	Light-curing glue	ADHX0001	PAA
F2	2	PE Bag	ADCX0011	PE
F3	1	2# Carton,,5pcs/pack, 20 packs/box	ADEX0001	Corrugated Paper

Customer Approval		LePure Biotech Co., Ltd.		
Company:		Drawing No. CR0500-1E		
Signature:		Title: SU Freezing Bags		
Date:		Designed By: Jiajia Guo	Date: 2023.03.17	Checked By: ShualkeLi
			Date: 2023.03.17	Edition: Rev03
Design Approval : Your signature on this drawing signifies your agreement to proceed with the manufacturing of the system(s) based on this design. By signing, you also assume full confidentiality and intellectual property. This design is the property of LePure Biotech Co., Ltd. and it contains confidential information that relates to our customer's process. We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property.				