

Pharma-Sat Plus™

Pharma-Sat Plus™ wipers are low endotoxin, ISO Class 5, non-woven, 45% polyester/55% cellulose blended wipers that are available in sterile and non-sterile versions. Our economical, non-shedding, pre-saturated wipers can be used throughout your clean-room for cleaning, disinfecting, spill control and wipe down of products. Pharma-Sat Plus™ wipers are pre-saturated with 70% Isopropanol and 30% USP WFI.



Thomas No. 20A00H630

Mfr. No. PSPW5VS797

PHARMA-SAT PLUS™ 70%
IPA Sterile Polyester/Cellulose
Wipers

Available in 9" x 9" sterile and
9" x 11" sterile.

SPECIFICATIONS

Size:	PSPW5VS797 (9" x 9")
Saturation:	70% IPA / 30% USP WFI
Composition:	45% polyester / 55% cellulose
Sterilization Method:	Gamma Irradiation
ISO Class:	ISO Class 5 and higher

Disclaimer: Pricing and offers subject to change without notice.