

# Driving science together

Boost lab productivity and save 40%\*

## Introduce automation into your workflow with Thermo Scientific™ sample storage starter packs

### Engineered with the value of your samples in mind

- Choose from five different starter pack options
- Eliminate the need to reformat non-2D cryotube samples
- Secure sample tracking with 2D, 1D, and human readable codes
- Cap/decap a full rack of tubes in less than 1 minute
- Scan 96 tubes in as little as 1 second with advanced and premium starter packs
- Help ensure quality with ISO-7 class cleanroom-manufactured tubes with sterility assurance





### Basic starter pack

Description	Quantity	Promotional price*	Cat. No.	Thomas No.
8-Channel Screw Cap Decapper	1 each	\$2,477 (Save \$1,651)	HHD-WP1D-2	23A00D886
0.5 mL Matrix ScrewTop Tri-coded Tubes	1 case			
1.0 mL Matrix ScrewTop Tri-coded Tubes	1 case			



### Advanced starter pack

Description	Quantity	Promotional price*	Cat. No.	Thomas No.
8-Channel Screw Cap Decapper	1 each	\$3,904 (Save \$2,602)	HHD-ST-WP1D-2	23A00D887
VisionMate ST Single-Tube 2D Barcode Reader	1 each			
0.5 mL Matrix ScrewTop Tri-coded Tubes	1 case			
1.0 mL Matrix ScrewTop Tri-coded Tubes	1 case			



### Premium starter pack

Description	Quantity	Promotional price*	Cat. No.	Thomas No.
8-Channel Screw Cap Decapper	1 each	\$11,109 (Save \$7,406)	HHD-HSX-WP1D-2	23A00D888
VisionMate HSX High Speed Barcode Reader	1 each			
0.5 mL Matrix ScrewTop Tri-coded Tubes	1 case			
1.0 mL Matrix ScrewTop Tri-coded Tubes	1 case			



### Premium 0.5 mL starter pack

Description	Quantity	Promotional price*	Cat. No.	Thomas No.
8-Channel Screw Cap Decapper	1 each	\$12,560 (Save \$8,372)	HHD-HSX-WP1D-505	23A00D889
VisionMate HSX High Speed Barcode Reader	1 each			
0.5 mL Matrix ScrewTop Tri-coded Tubes	5 cases			



### Premium 1.0 mL starter pack

Description	Quantity	Promotional price*	Cat. No.	Thomas No.
8-Channel Screw Cap Decapper	1 each	\$12,560 (Save \$8,372)	HHD-HSX-WP1D-510	23A00D890
VisionMate HSX High Speed Barcode Reader	1 each			
1.0 mL Matrix ScrewTop Tri-coded Tubes	5 cases			

\* This promotion is available only to life science professionals in the US and Canada who are 21 years or older. Discount will apply to qualifying orders received by Thermo Fisher Scientific no later than December 31, 2023, or until promotional supplies are depleted, whichever comes first. Customers can use the discount code an unlimited number of times. Discount applies to list price in effect at the time order is received by Thermo Fisher Scientific. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply.

**This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific use or application.** © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **COL119803 0323**

thermo scientific