

Universal Fit in Most Repeating Pipettor Brands

Roxy M™ Universal Fit Repeating Pipettor Tips



- Fit most manual repeating pipettors including Roxy M™, Eppendorf® Multipipette® 4780, Brand HandyStep® S, Fisherbrand® Repeating Pipettor, VWR® Repeating Pipettor, Gilson Distriman®, and Scilogex® StepMate Repeater Pipettor
- Internal plunger prevents liquid from contacting any part of the pipettor, eliminating the risk of contamination
- ISO 9001:2000 facility; Each pack of tips is lot stamped, comes with a Certificate of Conformity
- Sterile tips are individually packaged and bioclean® certified to be pyrogen-free, DNA-free, RNase-free, and ATP-free

TIPS	NO	1	2
0.05 ml		1 µl	2 µl
0.5 ml		10	20
1.25 ml		25	50
2.5 ml		50	100
5.0 ml		100	200
12.5 ml		250	500
25 ml		500	1000
50 ml		1000	2000
STEPS			23

Convenient chart of Dispensing Volumes per Roxy M tip size is provided on the pipettor.



Click Here
For a **FREE**
Sample Pack
of Roxy M™
Tips*

*FREE Sample Pack of Roxy M™ Tips

Includes:

1.25ml non-sterile, catalog no. 1209Q63

2.5ml non-sterile, catalog no. 1209Q64

5.0ml non-sterile, catalog no. 1209Q65



Adapter for use with 25 and 50ml tips

Catalog No.	Description
1209Q60	Roxy M™ Repeating Pipettor with 25/50ml tip adapter
1209Q61	Tips, 0.05ml, non-sterile
1209Q62	Tips, 0.5ml, non-sterile
1209Q63	Tips, 1.25ml, non-sterile†
1209Q64	Tips, 2.5ml, non-sterile†
1209Q65	Tips, 5.0ml, non-sterile†
1209Q66	Tips, 12.5ml, non-sterile
1209Q67	Tips, 25ml, non-sterile, requires adapter
1209Q68	Tips, 50ml, non-sterile, requires adapter
1209Q69	Tips, 0.5ml, sterile

†Included in sample pack

Catalog No.	Description
1209Q70	Tips, 1.25ml, sterile
1209Q71	Tips, 2.5ml, sterile
1209Q72	Tips, 5.0ml, sterile
1209Q73	Tips, 12.5ml, sterile
1209Q74	Tips, 25ml, sterile, requires adapter
1209Q75	Tips, 50ml, sterile, requires adapter
1209Q76	Tips, 0.05ml, sterile
1209Q77	Assortment Tip Pack - 5ea of 0.05ml-25ml non-sterile & 5ea 2.5ml & 5.0ml sterile
1209Q78	Adapters for 25ml and 50ml tips