

DTNB Polyclonal Antibody

Catalog number: 12045-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DTNB; DTNB, dystrobrevin,
beta, DTN B, Dystrobrevin beta,
Beta dystrobrevin



mouse skeletal muscle tissue were
subjected to SDS PAGE followed by
western blot with 12045-1-AP(DTNB
antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated DTNB MW:	560aa, 64 kDa
Observed DTNB MW:	64kd
Positive WB detected in	Mouse skeletal muscle tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2675
GenBank accession number:	BC016655
Gene ID (NCBI):	1838
Full name:	Dystrobrevin, beta

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.