

HIST3H2A Polyclonal Antibody

Catalog number: 10445-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

HIST3H2A; MGC3165



human brain tissue were subjected to SDS PAGE followed by western blot with 10445-1-AP(HIST3H2A Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue
Calculated HIST3H2A MW:	14kd
Observed HIST3H2A MW:	14-18kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0374
GenBank accession number:	BC001193
Gene ID (NCBI):	92815
Full name:	Histone cluster 3, H2a

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.