


Product Selection chart

		Capacity										Temp Accuracy	RS Interface	Analog Output	Temperature Alarm	Temperature Safety Device	Fan	Inner Glass Door	Controller	Power Outlet in Chamber
		20 liter = 0.7 cu ft	28 liter = 1 cu ft	53 liter = 1.9 cu ft	115 liter = 4.1 cu ft	150 liter = 5.4 cu ft	210 liter = 7.5 cu ft	240 liter = 8.5 cu ft	400 liter = 14.3 cu ft	700 liter = 24.7 cu ft	720 liter = 25.4 cu ft									
E	Oven with mechanical control																		Mechanical	
ED	Drying Oven w. gravity convection													V	3.1				Timer	LED
FD	Drying Oven w. forced convection											422	T	A	2				Timer/Ramp	LED
FED	Heating chamber											422	T	A	2	☉			Timer/Ramp	LED
FP	Temperature test chamber											422	T		2	☉			TF/Ramp	LED
FDL	Safety drying oven											422		A/V	2	☉			TF/Ramp	LED
MDL	Safety drying oven											422		A/V	2	☉			Color Graphic program	LED
VD	Vacuum drying oven											422		V	2				TF/program	LED
VDL	Safety vacuum drying oven											422		V	2				TF/program	LED
B	Incubator with mechanical control														1				Mechanical	
BD	Incubator											422	T	V	3.1				LED	T
BF	Incubator											422	T		3.1	☉			LED	TF/Ramp
KB	Refrigerated incubator											422	T	V/A	3.1	☉			LED	TF/program
CB	CO <sub>2</sub> incubator											422	T/CO <sub>2</sub>	V/A	3.1	☉			Color Graphic program	LED
C	CO <sub>2</sub> incubator												T/CO <sub>2</sub>	-	3.1	☉			LED	TF
KBW	Plant growth chamber											422	T	V/A	3.1	☉			LED	TF/program
KBWF	Plant growth chamber											422	T/H	V/A	3.1				Color Graphic program	LED
KBF	Constant climate chamber											422	T/H	V/A	3.3				Color Graphic program	LED
KBF P	Photostability test chamber											422		V/A	3.1				Color Graphic program	LED
KBF LQC	Constant climate chamber with LQC											422		V/A	3.1				Color Graphic program	LED
M	Temperature test chamber											422	T		2	☉			Color Graphic program	LED
MK	Environmental simulation chamber											422		V/A	2	☉			Color Graphic program	LED



A Audible  
 V Visual  
 H Humidity  
 E Ethernet  
 T Temperature  
 TF Timer functions  
 P Programmable

✓ Standard  
 ☉ Fan outside inner chamber  
 ☉ Speed control for fan  
 ■ Option