

SPINK1 Polyclonal Antibody

Catalog number: 13477-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SPINK1; PCTT, PSTI, SPINK1,
Spink3, TATI

Applications

Tested applications: ELISA
Species specificity: Human; other species not tested.
Calculated SPINK1 MW: 79aa,9 kDa
Observed SPINK1 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag4279
GenBank accession number: BC025790
Gene ID (NCBI): 6690
Full name: Serine peptidase inhibitor, Kazal type 1

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.