

## LGI1 Polyclonal Antibody

Catalog number: 12483-1-AP


Size: 33 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

LGI1; EPITEMPIN, Epitempin 1,  
EPT, ETL1, IB1099, LGI1



human brain tissue were subjected to SDS PAGE followed by western blot with 12483-1-AP(LGI1 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated LGI1 MW:	557aa,64 kDa
Observed LGI1 MW:	64 kDa,58 kDa
Positive WB detected in	Human brain tissue, HeLa cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3184
GenBank accession number:	BC022500
Gene ID (NCBI):	9211
Full name:	Leucine-rich, glioma inactivated 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.