

## PPME1 Polyclonal Antibody

Catalog number: 14435-1-AP


Size: 27 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PPME1; FLJ22226; PME-1



A549 cells were subjected to SDS PAGE followed by western blot with 14435-1-AP (PPME1 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	A549 cells, HeLa cells, human brain tissue
Calculated PPME1 MW:	42kd
Observed PPME1 MW:	44kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag5817
GenBank accession number:	BC050705
Gene ID (NCBI):	51400
Full name:	Protein phosphatase methylesterase 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.