

PKP4 Polyclonal Antibody

Catalog number: 14062-1-AP


Size: 34 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PKP4; FLJ31261; FLJ42243;
p0071



A431 cells were subjected to SDS PAGE followed by western blot with 14062-1-AP (PKP4 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	A431 cells, HepG2 cells
Calculated PKP4 MW:	132kd
Observed PKP4 MW:	125-140kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5113
GenBank accession number:	BC050308
Gene ID (NCBI):	8502
Full name:	Plakophilin 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.