

## C17orf81 Polyclonal Antibody

Catalog number: 10162-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

C17orf81; C17orf81, DERP6,

ELP5, HSPC002, MST071,

MSTP071, S phase 2 protein, S-phase 2 protein

### Applications

Tested applications: ELISA

Species specificity: Human; other species not tested.

Calculated C17orf81 MW: 34.8 kDa

Observed C17orf81 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen: Ag0194

GenBank accession number: BC002762

Gene ID (NCBI): 23587

Full name: Chromosome 17 open reading frame 81

### Product information

Purification method: Antigen affinity purification

Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.