

MAPK7 Polyclonal Antibody

Catalog number: 11164-1-AP


Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MAPK7; Big MAP kinase 1,
BMK 1, BMK1, ERK 5, ERK4,
ERK5, MAP kinase 7, MAPK 7,
MAPK7, PRKM7



SKOV-3 cells were subjected to SDS PAGE followed by western blot with 11164-1-AP (MAPK7 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated MAPK7 MW:	89 kDa
Observed MAPK7 MW:	70 kDa
Positive WB detected in	SKOV-3 cells, HEK-293 cells, mouse heart tissue, mouse testis tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1648
GenBank accession number:	BC009963
Gene ID (NCBI):	5598
Full name:	Mitogen-activated protein kinase 7

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.