


- Sizes:** 21"/24"/26"
- Thomas No.:** 20A00J384 (XC2100-B5), 20A00J385 (XC2100-W5), 20A00J386 (XC2400-B5), 20A00J387 (XC2400-W5), 20A00J388 (XC2100-W1), 20A00J389 (XC2400-B1), 20A00J390 (XC2400-W1), 20A00J391 (XC2600-W1).
- Packaging:** 100pcs/bag,1000pcs/case
- Colors:** White/blue
- Material:** Non-woven spunbonded polypropylene, 10gsm
- Elastic:** Latex-free
- Specifications:**

### XC2000 Material

#### Physical Properties

Test Item	Test Standard	Result
Weight (g/m <sup>2</sup> )	ASTM D3776/D-3776M-20 (Option C)	11.7
Tensile Strength Length, lbf Width, lbf	ASTM D 5034-09, Grab Method	5.0 4.2



Item	21"	24"	26"	Tolerance	Comment
A	21cm	24cm	26cm	+/-2cm	Elastic relaxed
B	53cm	61cm	66cm	+/-2cm	Stretched length