

PRMT8 Polyclonal Antibody

Catalog number: 12557-1-AP


Size: 27 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PRMT8; HRMT1L3; HRMT1L4



MCF7 cells were subjected to SDS PAGE followed by western blot with 12557-1-AP (PRMT8 Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human brain using 12557-1-AP (PRMT8 Antibody) at Dilution 1:50 (under 40x lens)

Applications

| | |
|---------------------------|--|
| Tested applications: | ELISA, WB, IHC |
| Species specificity: | Human, mouse, rat; other species not tested. |
| Positive WB detected in: | MCF7 cells, Jurkat cells, SH-SY5Y cells |
| Calculated PRMT8 MW: | 394aa, 45kd |
| Observed PRMT8 MW: | 45kd |
| Positive IHC detected in: | Human brain |
| Recommended dilution: | WB: 1:1000-1:10000 IHC: 1:20-1:200 |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

| | |
|---------------------------|--------------------------------------|
| Immunogen: | Ag3166 |
| GenBank accession number: | BC022458 |
| Gene ID (NCBI): | 56341 |
| Full name: | Protein arginine methyltransferase 8 |

Product information

| | |
|----------------------|--|
| Purification method: | Antigen affinity purification |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |