

RBM14 Polyclonal Antibody

Catalog number: 10196-1-AP

Size: 50 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

RBM14; COAA;

DKFZp779J0927; MGC15912;

MGC31756; PSP2; SIP; SYTIP1;

TMEM137



Jurkat cells were subjected to SDS PAGE followed by western blot with 10196-1-AP (RBM14 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Jurkat cells
Calculated RBM14 MW:	70kd
Observed RBM14 MW:	70-75kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0254
GenBank accession number:	BC000488
Gene ID (NCBI):	10432
Full name:	RNA binding motif protein 14

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.