



Cole-Parmer® VPS-200 Series Air Vacuum Pumps

- Compact size and integral handle make pumps easily transportable
- Economical vacuum pumps serve as a pressure or vacuum source with an excellent price to performance ratio
- Low-friction piston runs cool and quiet for a long service life with consistent performance
- 1-year warranty



Vacuum/Pressure Pump
79202-00 | 79202-05

Cole-Parmer VPS-200 Series Air Vacuum Pumps (including Stations)

This portable vacuum pump is designed for portability and versatility with a tough, lightweight plastic housing and integral carrying handle. The long-wearing motor delivers performance characteristics suitable for a wide range of applications.

What's included: 1.8 m (6 ft) cord with US standard plug on 115 VAC models. IEC socket with 1.8 m (6 ft) EU and UK mains lead with plug on 230 VAC models.

NOTE: Do not pump liquids or immerse pump. Avoid organic, basic and acidic vapors.



Vacuum Pump Station
79202-30 | 79202-35

Specifications

Voltage	Maximum vacuum	Maximum pressure	Port connections	Motor	Max temperature	Humidity range	Item number
115 VAC	541 (±0.2) mm Hg 21.3 (±0.01) in Hg	33 (±0.01) psig* 22.8X10 ⁴ Pa	1/4" hose barb	1/11.5 hp (65w) shaded pole	40 °C (104 °F)	10% to 90% non-condensing	79202-00
230 VAC							79202-05
115 VAC	541 (±0.2) mm Hg 21.3 (±0.1) in Hg	33 (±0.01) psig* 22.8X10 ⁴ Pa	1/4" hose barb	1/11.5 hp (65w) shaded pole	40 °C (104 °F)	10% to 90% non-condensing	79202-30
230 VAC							79202-35

*Intermittent duty: limit output to 33 psig – Continuous duty: limit output to 10 psig.

Ordering Information

Description	Dimensions (L x W x H)	Weight	Electrical	Item number
VPS-200 Air Diaphragm Vacuum/Pressure Pump, 0.37 cfm	18.4 x 10.9 x 21.5 cm (7.2 x 4.3 x 8.5 in)	1.6 kg (3.5 lb)	115 VAC, 60 Hz	79202-00
VPS-200 Air Diaphragm Vacuum/Pressure Pump, 0.31 cfm			230 VAC, 50 Hz	79202-05
VPS-200 Vacuum Pump Station with Pressure Gauge and Control, 0.37 cfm	23.05 x 17.87 x 21.65 cm (9.07 x 7.03 x 8.52 in)	2.8 kg (6.17 lb)	115 VAC, 60 Hz	79202-30
VPS-200 Vacuum Pump Station with Pressure Gauge and Control, 0.31 cfm			230 VAC, 50 Hz	79202-35

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

xxxx

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480.272279
All other countries: +1.847.549.7600

Cole-Parmer®