

# KapSafe-mini™


## Bench-Top Automated Recapper



**Protect Against Repetitive -Motion  
Injuries Associated with Manual  
Recapping**



- Bench-top system for the recapping of single racks
- Recaps original vacuum collection samples; 13 & 16 mm diameter, 75-100 mm high & the LGP 12x75mm transfer tubes
- Uses unique LGP caps for multiple tube sizes with capacity of up to 1,000 caps
- Handles from 5 position to 32 position single racks, inc. LGP generic 10 pos. rack
- Recaps up to 1,200 tubes per hour
- Minimizes evaporation and sample contamination for stored sample tubes
- Pneumatics-free, small footprint (22" w. x 14" d. x 26" h.), 120VAC & 220VAC



Extensions  
for longer  
racks are  
available

**KapSafe-mini™ IS THE ANSWER!**

- STOP Manually Recapping Tubes!
- Protect Against Repetitive Motion Injuries
- Protect Against Tech Boredom!
- Improve Laboratory Productivity