

TSPYL1 Polyclonal Antibody

Catalog number: 13932-1-AP


Size: 39 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TSPYL1; TSPYL



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13932-1-AP(TSPYL1 Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse; other species not tested.
Positive WB detected in:	Mouse lung tissue, mouse heart tissue, mouse kidney tissue
Caculated TSPYL1 MW:	49kd
Observed TSPYL1 MW:	49 kDa
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4961
GenBank accession number:	BC048969
Gene ID (NCBI):	7259
Full name:	TSPY-like 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.