

EISCO LAB STARTER KIT

ESSENTIAL METALWARE AND UTENSILS TO OUTFIT YOUR LAB BENCH IN ONE CONVENIENT PACKAGE



Complete, 13-piece kit with chemically-resistant, powder-coated alloy racks and supports includes:

- Retort Base, 8"x5", with stainless steel rod, 23.6"L
- Universal Clamp/Bosshead, 10cm max opening
- Ring Clamp/Bosshead, 7.6cm I.D.
- Double Burette Clamp, holds up to 100mL burettes
- Tripod Stand, 19cm tall x 7.5cm I.D.
- Wire Mesh Square, 12cm sides
- 17-place Test Tube/Utensil rack, 12-25mm dia. holes
- Test Tube Brush
- Test Tube Tongs
- Flask Tongs
- Crucible Tongs
- Lab Spoon/Spatula

Special Sale Price!

NOW

Original Price ~~\$126.20~~

\$89⁹⁹

Eisco

Lab Starter Kit

[Thomas No. 1182Z85](#)

Special pricing promotion valid
April 1 through July 31, 2024

Disclaimer: Pricing and offers subject to change without notice. 3/28/24