

## Specification Sheet

**PRODUCT: Purified Water, USP**

CATALOG NO.:	521009P	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	12 Months from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C	CONTAINER SIZE:	200L

**Test Results**

Test	Specification	Test Method	Result
Conductivity	Meets Requirements ( $\leq 1.3 \mu\text{S}/\text{cm}$ @ 25°C)	USP <645> Bulk Water	TBD
Total Organic Carbon	$\leq 500$ ppb	USP <643> Bulk Water	TBD
Heterotrophic Plate Count	$\leq 100$ CFU/mL	USP <1231>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

**Manufacturing Site**

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538