

For Research Use Only



# TTC26 POLYCLONAL ANTIBODY

Catalog Number: 25083-1-AP

## Basic Information

Catalog Number:  
25083-1-AP

Size:  
30 µg/150 µl

Source:  
Rabbit

Isotype:  
IgG

Purification Method:  
Antigen affinity purification

Immunogen Catalog Number:  
AG18983

GenBank Accession Number:  
BC126331

GeneID (NCBI):  
79989

Full Name:  
tetra(ricopeptide repeat domain 26

Calculated MW:  
554aa, 64 kDa

Observed MW:  
60-64 kDa

Recommended Dilutions:  
WB 1:200-1:1000

## Applications

Tested Applications:  
ELISA, WB

Species Specificity:  
Human, Mouse

Positive Controls:

WB : mouse testis tissue; HEK-293 cells, NIH/3T3 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 25083-1-AP( TTC26 Antibody) at dilution of 1:300