



## AATCC Mineral Salts Agar

M232

AATCC Mineral Salts Agar is recommended for evaluation of fungicides used on textiles and to study the resistance of textiles to mildew and rot.

### Composition\*\*

Ingredients	Gms / Litre
Ammonium nitrate	3.000
Monopotassium phosphate	2.500
Dipotassium phosphate	2.000
Magnesium sulphate	0.200
Ferrous sulphate	0.100
Agar	20.000
Final pH ( at 25°C)	5.6±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 27.8 grams in 1000 ml distilled water. Add 7.5 grams of dextrose if the medium is to be used for testing with *Aspergillus brasiliensis*. Heat to boiling to dissolve the medium completely. Dispense in 7 ml amounts in test tubes or 40 ml amounts in bottles or flasks. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

AATCC Mineral Salts Agar is used as per the procedure described by American Association of Textile Chemists and Colourists (1). These procedures are employed for testing a) Initial mildew resistance of textiles, b) Perseverance of mildew resistance, c) Fungicidal potency of textile fungicides testable, Standard Minimum protective concentration, peccnanence indices, including resistance to leaching, volatalization and weathering. Cultures used are *Chaetomium globosum* ATCC 16790, *Aspergillus brasiliensis* ATCC 16404 (2).

Ammonium nitrate acts as a nitrogen source. Dipotassium and monopotassium phosphate provides buffering to the medium. Magnesium sulphate and ferrous sulphate are sources of ions that simulate metabolism.

### Quality Control

#### Appearance

Cream to beige homogeneous free flowing powder

#### Gelling

Firm, comparable with 2.0% Agar gel.

#### Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel.

#### Reaction

Reaction of 2.78% w/v aqueous solution at 25°C. pH : 5.6±0.2

#### pH

5.40-5.80

#### Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

#### Cultural Response

Organism	Inoculum (CFU)
----------	----------------

#### Cultural Response

* <i>Aspergillus brasiliensis</i> ATCC 16404	50-100
---	--------

*Chaetomium globosum* 50-100  
ATCC 16790

### Storage and Shelf Life

Store below 30°C and the prepared medium at 2-8°C. Use before expiry date on the label.

### Reference

1. Technical Manual of the American Association of Textile Chemists and Colourists, 1959, Part III, AATCC Test Methods, 82, Lowell, Mass
2. Catalogue of Bacteria and Bacteriophages, 1992, 18th Ed., American Type Culture Collection, Rockville, MD

Revision : 1 / 2011



#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.