

CELL STRAINERS PROMOTION

BUY 2 GET 1 FREE!



Cell Strainers can be used for many applications, including isolating cells and removing clumps or debris

- Unique extended lip on housing enables aseptic handling with forceps
- Elongated gripping tab on housing for easy handling
- Hydrophilic treatment increases flow and reduces the effects of surface tension
- 3 Mesh size options, with color-coding for easy identification
- Universally fits most brands of 50mL conical tubes
- Exclusive Pestle improves manipulation of material in the strainer

	Part No.	Product Description	Qty/Case
Cell Strainers — 15mL			
NEW	230471	15mL Cell Strainer, 40µm, Blue	50
NEW	230473	15mL Cell Strainer, 70µm, White	50
NEW	230475	15mL Cell Strainer, 100µm, Yellow	50
Cell Strainers — 50mL			
	229481	50mL Cell Strainer, 40µm, Blue	50
	229482	50mL Cell Strainer, 40µm, Blue, Bulk, Pk	50
	229483	50mL Cell Strainer, 70µm, White	50
	229484	50mL Cell Strainer, 70µm, White, Bulk, Pk	50
	229485	50mL Cell Strainer, 100µm, Yellow	50
	229486	50mL Cell Strainer, 100µm, Yellow, Bulk, Pk	50
Cell Strainers — 500mL			
NEW	230491	500mL Cell Strainer, 40µm, Blue	50
NEW	230493	500mL Cell Strainer, 70µm, White	50
NEW	230495	500mL Cell Strainer, 100µm, Yellow	50
Cell Strainers — Accessories			
	229480	Pestle, Green	50

Buy 2 of same product to get 1 Free of your choice

How To Redeem:

OPTION 1:



Submit online at:
celltreat.us/submitpromo

PROMO CODE:

B2G1-STRAINERS-TSI26-Q3

Offer valid in the U.S. only and void where prohibited by government, law or company policy. Offer valid for end-user customers only for orders placed between July 1, 2026 and September 30, 2026. Limit 4 free cases per lab. Allow 4 - 6 weeks for delivery of product. Offer does not apply to previous purchases and cannot be combined with any other CELLTREAT/VistaLab discounts or promotional offers. Submission form must be received by October 15, 2026.



www.ThomasSci.com
CustomerService@ThomasSci.com

OPTION 2:

Email form & proof of purchase to
promotions@celltreat.com

Name: _____

Company: _____

Shipping Address: _____

Lab/Room #: _____

City: _____ State: _____ Zip: _____

Email: _____

FREE CASE SELECTION - Eligible part numbers above, limit 4 FREE cases per lab

Case 1		Case 3	
Case 2		Case 4	