


Weighing Supplies

A variety of styles to suit your weighing requirements

Aluminum Moisture Pans

Ideal for loss and drying determinations

- Withstand temperatures up to 677°C (1,250°F)
- Can be used as evaporating dishes or dust covers
- Offered with a flat bottom, slightly crimped side walls and a small tab handle, or with smooth walls and no tab handle



9885D01

9885D19

Thomas No.	Size (mm)	Diam. (In.)	Approx. Cap.	Depth (mm/in.)	Weight	Pack	Case
Aluminum Weigh Dishes with Crimped Side Walls and Tab Handle							
9885D01	43	1.75	20 mL	7 mm / 5/16 in.	1.0 g	100	1000
9885D07	57	2.25	50 mL	16 mm / 5/8 in.	1.5 g	100	1000
9885D13	70	2.75	70 mL	17 mm / 11/16 in.	2.0 g	100	1000
Aluminum Weigh Dishes with Smooth Walls and No Tab Handle							
9885D19	57	2.25	50 mL	16 mm / 5/8 in.	1.5 g	100	1000
9885D28	70	2.75	70 mL	17 mm / 11/16 in.	2.0 g	100	1000

Pre-Cut Aluminum Squares *Eliminate time spent cutting aluminum foil*

- Multiple uses: Cap flasks, beakers, bottles and other labware
- Act as a barrier or insulator
- Easy to separate: Squares are interleaved with paper
- High resistance to fats, petroleum, greases and organic solutions non-toxic, non-magnetic, non-sparking and corrosion resistant; good electrical and thermal conductivity



Continued on reverse side

Thomas No.	Size (In.)	Thickness (In.)	No. of Squares*
1086J24	4 x 4	0.001	2200
1086J26	6 x 6	0.001	1100
1086J28	8 x 8	0.0015	430
1086J30	6 x 6	0.0015	192
1203N63	12 x 12	0.0015	1000

*Furnished in 5 lb. packages; Count is approximate

Weighing Dishes

- Available in Hexagonal and Square
- Flat bottom
- Perfect for weighing solids and liquids
- Made from flexible polystyrene
- Provides excellent protection against spills
- Prevents hazards and sample loss



Thomas No.	Size	Dimensions (In.)	Pack	Case
Square				
9885D37	39 mm	1-1/2 x 1-1/2 x 3/8	500	6000
9885D40	77 mm	3 x 3 x 1	500	2000
9885D43	127 mm	5 x 5 x 1	500	1000



Thomas No.	Dimensions (In.)	Box	Case
Hexagonal			
9885D55	1-1/2 x 1	500	4000
9885D58	2-1/2 x 2	500	3000
9885D60	4 x 3	500	2000
9885D63	5 x 3-1/2	500	2000




Weighing Papers

- Economical
- Smooth, strong and non-absorbent with glossy surface
- Especially suited for weighing powders

Thomas No.	Dimensions (In.)	Box
9885G40	3 x 3	500
9885G50	4 x 4	500
9885G60	6 x 6	500

Weighing Canoes



- Aerospace-grade, anti-static plastic
- Natural opaque white color
- Economical, stable, easy-to-handle design avoids "squeezing" of square weigh boats
- Plastic provides reduced adhesion without oily coatings

Thomas No.	Dimensions (In.)	Pack	Case
9885E07	2-1/4 x 1 x 1/4	250	3000
9885E09	5 x 3-1/4 x 1	250	1500
9885E11	7 x 4-1/4 x 1-1/8	250	1000