



## CERTIFICATE OF ANALYSIS

### Human Albumin 25% Solution

(Catalog Number 800-120)

Lot Number: H122031

Date of Manufacture: Jun2022

Product Expiry: Nov2023

Origin: United States

Storage Temperature:  $\leq -10^{\circ}\text{C}$

Test	Methodology	Specification	Analysis
<b>Microbiological Testing</b>			
Sterility	USP<71>		
Bacteria		No Growth	No Growth
Fungi		No Growth	No Growth
<b>Viral Testing</b>			
	21 CFR 610.40		
HBsAg		Non-Reactive	Non-Reactive
Anti-HCV		Non-Reactive	Non-Reactive
Anti-HIV-1/HIV-2		Non-Reactive	Non-Reactive
Syphilis		Negative	Negative
HIV-1 NAT		Non-Reactive	Non-Reactive
HBV-NAT		Non-Reactive	Non-Reactive
HCV-NAT		Non-Reactive	Non-Reactive
<b>Physical Testing</b>			
Purity		Test and Report	100%
Protein		Test and Report	23.0 g/dL
Endotoxin		Test and Report	1.32 EU/mL
pH		Test and Report	6.99

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use. Results shown were obtained by carefully performed methods believed to be reliable. However, since some results may vary for specific tests depending upon methodology and other variables, it is suggested that tests for which results are particularly important be repeated by the user of this product.

All human plasma is collected from stringently screened donors at FDA-licensed collection centers in the United States.

Jaqueline Wilson  
Name

QA Specialist  
Title

20 JUL 2022  
Date