

APOM Polyclonal Antibody

Catalog number: 12817-1-AP

Size: 28 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:


APOM; Apo M, apolipoprotein

M, APOM, G3a, HSPC336,

NG20, Protein G3a



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using 12817-1-AP(APOM Antibody) at dilution of 1:50 (under 40x lens)



human plasma tissue were subjected to SDS PAGE followed by western blot with 12817-1-AP(APOM Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated APOM MW:	21 kDa
Observed APOM MW:	25 kDa
Positive WB detected in	Human plasma tissue
Positive IHC detected in	Human hepatocirrhosis tissue
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3137
GenBank accession number:	BC020683
Gene ID (NCBI):	55937
Full name:	Apolipoprotein M

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.