

GEB Universal tips for TWC thin wallcrown And for improved flexibility.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (3.8*3.5)mm	Packaging Spec.
FT0100-B-N	100ul	Natural			Filtered	1000tips/bag,
FT0100-B-N-LB	100ul	Natural	Low Retention		Filtered	10bags/case,
FT0100-B-NS	100ul	Natural		Sterilized	Filtered	10000tips/case
FT0100-B-NS-LB	100ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 480*345*280mm.			Gross Weight: 4.5kg.			
FT0100-R-NS	100ul	Natural		Sterilized	Filtered	96tips/rack,
FT0100-R-NS-LB	100ul	Natural	Low Retention	Sterilized	Filtered	10racks/pack, 5packs/case, 4800tips/case. 50racks/case.
Case dimension L*W*H: 500*265*450mm.			Gross Weight: 7.8kg.			

- ✧ 100μL with graduation marks at 10ul,50ul
- ✧ 96 Tip/Hanging Tip Rack format / Bulk
- ✧ Polypropylene
- ✧ Presterile/Nonsterile
- ✧ Filtered (With aerosol-resistant filters) /NonFiltered
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Universal tips workstations.

Drawing: 100ul

