

UBTD2 Polyclonal Antibody

Catalog number: 13172-1-AP


Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

UBTD2; DCUBP, DC-UbP,
ubiquitin domain containing 2,
Ubiquitin like protein SB72,
UBTD2



Immunohistochemical of paraffin-embedded human breast cancer using 13172-1-AP (UBTD2 antibody) at dilution of 1:100 (under 10x lens)



human liver tissue were subjected to SDS PAGE followed by western blot with 13172-1-AP (UBTD2 antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated UBTD2 MW:	234aa, 26kd
Observed UBTD2 MW:	26kd
Positive WB detected in:	Human liver tissue
Positive IHC detected in:	Human breast cancer tissue
Recommended dilution:	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3827
GenBank accession number:	BC019910
Gene ID (NCBI):	92181
Full name:	Ubiquitin domain containing 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.