

Solvent dispenser

Ideal for a wide ranges of applications in the electronic industry as well as for chemistry and medical laboratories. Made of anti-static HDPE plastic with Teflon tube inside, for a great chemical resistance.

The main advantages are:

- » dispenses clean
- » prevents evaporation and contamination
- » suitable for most liquids like solvents, degreasers, oils
- » easy use
- » reduces waste
- » prevents fumes, accident and spillage
- » ESD safe
- » surface resistivity 10¹⁰ Ohm
- » no charge generation, no fire or explosion ignition risk

No stopper or cap to remove. Just touch the dispensing dish with a brush, cloth or other applicator and solvent is brought to the surface the instant it is needed. Stainless steel valves seal contents in bottle until used, preventing spills, dripping and escaping fumes from endangering the technician. Fluid pumped into the dish does not drip back into the bottle, keeping stored liquids uncontaminated. Bottles are available in a variety of sizes.



TSD12



TSD12AS



TSD18



TSD18AS

Art. Nr.	Description
TSD12	120 ml safety solvent dispenser - pink
TSD12AS	120 ml safety solvent dispenser with anti-splash system - pink
TSD18	180 ml safety solvent dispenser - pink
TSD18AS	180 ml safety solvent dispenser with anti-splash system - pink