

MVK Polyclonal Antibody

Catalog number: 12228-1-AP


Size: 30 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

MVK; FLJ96772, LRBP,
mevalonate kinase, MK, MVK,
MVLK



mouse uterus tissue were subjected to SDS PAGE followed by western blot with 12228-1-AP(MVK antibody) at dilution of 1:500



Immunofluorescent analysis of PC-3 cells using 12228-1-AP(MVK Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Applications

Tested applications:	ELISA, WB, IF
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated MVK MW:	396aa,42 kDa
Observed MVK MW:	42 kDa
Positive WB detected in	Mouse uterus tissue, PC-3 cells
Positive IF detected in	PC-3 cells
Recommended dilution:	WB: 1:500-1:5000 IF: 1:10-1:100

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2865
GenBank accession number:	BC016140
Gene ID (NCBI):	4598
Full name:	Mevalonate kinase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.