

LY6K Polyclonal Antibody

Catalog number: 12026-1-AP

Size: 27 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:


LY6K; CO16, FLJ35226,

HSJ001348, Ly 6K, LY6K,

Lymphocyte antigen 6K



human testis tissue were subjected to SDS PAGE followed by western blot with 12026-1-AP(LY6K antibody) at dilution of 1:500



IP Result of anti-LY6K (IP:12026-1-AP, 4µg; Detection:12026-1-AP 1:300) with P C-3 cells lysate 3200ug.

Applications

Tested applications:	ELISA, WB, IP
Species specificity:	Human; other species not tested.
Calculated LY6K MW:	179aa,20 kDa
Observed LY6K MW:	28 kDa
Positive WB detected in	Human testis tissue
Positive IP detected in	PC-3 cells
Recommended dilution:	WB: 1:200-1:2000 IP: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2656
GenBank accession number:	BC001291
Gene ID (NCBI):	54742
Full name:	Lymphocyte antigen 6 complex, locus K

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.