

ROTORS WITH DIFFERENT SPEED/G-FORCE DEPENDING ON CENTRIFUGE

| Model: | Z306 | Z326 | Z326-K | Z366 | Z366-K |
|--------------------------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Thomas Part No Rotor Part No. | | | | | |
| 1159P76 Z326-100H-AC | 4,000 rpm/ 2,486 xg | 5,000 rpm/ 3,885 xg | 5,000 rpm/ 3,885 xg | | |
| 1159P84 Z326-2420 | 14,000 rpm/ 18,624 xg | 14,000 rpm/ 18,624 xg | 15,000 rpm/ 21,379 xg | 15,000 rpm/ 21,379 xg | 16,000 rpm/ 24,325 xg |
| 1159P85 Z326-2420H | 14,000 rpm/ 18,624 xg | 14,000 rpm/ 18,624 xg | 15,000 rpm/ 21,379 xg | 15,000 rpm/ 21,379 xg | 16,000 rpm/ 24,325 xg |
| 1159P87 Z326-3020H | 12,000 rpm/ 15,131 xg | 13,000 rpm/ 17,758 xg | 14,000 rpm/ 20,505 xg | | |
| 1159P89 Z326-4420 | | 15,000 rpm/ 21,379 xg | 15,000 rpm/ 21,379 xg | 15,000 rpm/ 21,379 xg | 16,000 rpm/ 24,325 xg |
| 1159P82 Z326-1220H | 18,000 rpm/ 23,542 xg | | | | |
| 1154K71 Z326-1205 | | 14,000 rpm/ 18,624 xg | 15,000 rpm/ 21,379 xg | | |
| 1159P83 Z326-2010 | | 12,000 rpm/ 15,775 xg | 12,000 rpm/ 15,775 xg | 12,000 rpm/ 15,775 xg | 12,000 rpm/ 15,775 xg |
| 1159P69 Z326-0650-HS | | 12,000 rpm/ 13,522 xg | 12,000 rpm/ 13,522 xg | 13,000 rpm/ 15,869 xg | 13,000 rpm/ 15,869 xg |
| 1159P79 Z326-1050 | | 7,500 rpm/ 8,174 xg | 7,500 rpm/ 8,174 xg | 9,000 rpm/ 11,771 xg | 10,500 rpm/ 16,022 xg |
| 1159P66 Z326-0485 | | | | | |
| 1159P71 Z326-0685 | | 11,000 rpm/ 15,555 xg | 13,000 rpm/ 21,726 xg | | |
| 1159P72 Z326-0685-25 | | 11,000 rpm/ 13,932 xg | 13,500 rpm/ 20,984 xg | | |
| 1159Q13 Z446-04500 | | | | | |

| Model: | Z446 | Z446-K | Z496-K | Z496-K-UC | Z32-HK |
|--------------------------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Thomas Part No Rotor Part No. | | | | | |
| 1159P76 Z326-100H-AC | | | | | 5,000 rpm/ 3,885 xg |
| 1159P84 Z326-2420 | 16,000 rpm/ 24,325 xg | 16,000 rpm/ 24,325 xg | | | 17,000 rpm/ 27,460 xg |
| 1159P85 Z326-2420H | 16,000 rpm/ 24,325 xg | 16,000 rpm/ 24,325 xg | | | 20,000 rpm/ 38,007 xg |
| 1159P87 Z326-3020H | | | | | 17,000 rpm/ 30,368 xg |
| 1159P89 Z326-4420 | 16,000 rpm/ 24,325 xg | 16,000 rpm/ 24,325 xg | | | 15,000 rpm/ 21,379 xg |
| 1159P82 Z326-1220H | | | | | |
| 1154K71 Z326-1205 | | | | | |
| 1159P83 Z326-2010 | 12,000 rpm/ 15,775 xg | 12,000 rpm/ 15,775 xg | | | 14,000 rpm/ 21,472 xg |
| 1159P69 Z326-0650-HS | 13,000 rpm/ 15,869 xg | 13,000 rpm/ 15,869 xg | | | 16,000 rpm/ 24,039 xg |
| 1159P79 Z326-1050 | 10,000 rpm/ 14,532 xg | 10,500 rpm/ 16,022 xg | 10,500 rpm/ 16,022 xg | 10,500 rpm/ 16,022 xg | 9,000 rpm/ 11,770 xg |
| 1159P66 Z326-0485 | 15,000 rpm/ 23,140 xg | 16,000 rpm/ 26,328 xg | | | 15,000 rpm/ 23,140 xg |
| 1159P71 Z326-0685 | 11,000 rpm/ 15,555 xg | 13,000 rpm/ 21,726 xg | | | |
| 1159P72 Z326-0685-25 | | | | | 13,000 rpm/ 21,140 xg |
| 1159Q13 Z446-04500 | 8,000 rpm/ 10,367 xg | 8,000 rpm/ 10,367 xg | 8,000 rpm/ 10,374 xg | 8,000 rpm/ 10,374 xg | |