

## Calicin Polyclonal Antibody

Catalog number: 13972-1-AP


Size: 31 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CCIN; calicin, CCIN



human testis tissue were subjected to SDS PAGE followed by western blot with 13972-1-AP(CCIN antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated Calicin MW:	67 kDa
Observed Calicin MW:	60-65kd
Positive WB detected in	Human testis tissue, rat testis tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5005
GenBank accession number:	BC019251
Gene ID (NCBI):	881
Full name:	Calicin

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.