

SIAH2 Polyclonal Antibody

Catalog number: 12651-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SIAH2; hSiah2, Seven in
absentia homolog 2, Siah 2,
SIAH2



mouse brain tissue were subjected to SDS PAGE followed by western blot with 12651-1-AP(SIAH2 Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated SIAH2 MW:	324aa, 35 kDa
Observed SIAH2 MW:	35-40 kDa
Positive WB detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3342
GenBank accession number:	BC013082
Gene ID (NCBI):	6478
Full name:	Seven in absentia homolog 2 (Drosophila)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.