



Item No.	Qty	Description	Inventory Code	Material
D1	1	5L SU Storage Bags	ACAX0012	ULDPE
B1	2	3/8" 5/8" (9.6*15.9mm) L-Flex Tubing, Length 500mm	DBBX0006	TPE
B2	1	1/4" 3/8" (6.4*9.6mm) Silicone Tubing (LePure), Length 200mm	DAAX0066	Silicone
C1	1	1/4" (6.4mm) Female Luer Thread Barb	EAXX0143	PP
C2	1	Needle-less sampler	EABX0040	PC
C3	1	Needle-less sampler cap	EABX0005	PC
C4	1	3/8" (9.6mm) Hose Barb Non-Valved In-Line Coupling Insert	EBAX0021	PC
C5	1	MPC Coupling Cap	EBAX0025	PP
C6	1	3/8" (9.6mm) Hose Barb Non-Valved In-Line Coupling Body	EBAX0009	PC
C7	1	MPC Coupling Plug	EBAX0028	PC
F1	2	2.5mm cable tie	ADAA0005	Nylon
F2	7	3.6mm cable tie	ADAA0004	Nylon
F3	1	Robert clip (small)	EDBA0026	PA
F4	2	Pipe clamp	EDBA0027	PP
F5	4	Silica gel end cap (medium clip)	ADBX0003	Silicone

<b>Customer Approval</b>		<b>LePure Biotech Co., Ltd.</b>				
Company:		Drawing No. <b>HP005-5AK</b>				
Signature:		Title: <b>SU Storage Bags</b>				
Date:		<b>Designed By:</b>	<b>Date:</b>	<b>Checked By:</b>	<b>Date:</b>	<b>Edition</b>
		Jiajia Guo	2023.03.12	Alan	2023.03.12	Rev01
<small>Design Approval the manufacturing of the system(s) based on this design. By signing, you also assume full</small>		<small>We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property</small>				